

Drager evita XL

- Adult Pediatric (Neo Optional Available)
- High End Premium ICU Ventilator
- PPV / assist, CMV / Assist, MMV+PS, APRV, BiPAP & CPAP
- Graphics, Loops, Trends & Auto Flow
- Lung Mechanics & Expiratory Parameters



Drager Evita 4

- High End Premium ICU Ventilator
- Adult Pediatric (Neo Optional Available)
- PPV / assist, CMV / Assist, MMV+PS, APRV, BiPAP & CPAP
- 10" color FTP Touch Screen
- Graphics, Loops, Trends & Auto Flow
- Lung Mechanics & Expiratory Parameters



Siemens Servo 300 / 300A

- Neonate To Adult Ventilator
- Volume & Pressure Control Modality
- Monitoring of All Exhaled Parameters
- Graphical Screen For Graphs And Loops
- 300AIS with auto Mode
- Compliance & Resistance Display



Drager - Babylog - 8000 Plus

- Dedicated Neonatal Ventilator
- Pressure Control with Volume Guarantee
- Monitoring of All Exhaled Parameters AT ET Level
- Accuracy of 0.1 mL
- Graphical Display for Graph
- Compliance & Resistance Display



Pesptronics - ESPRIT

- BiPAP (NIV) Mode
- Paediatric-adult Ventilator
- Blower driven Machine with
- Touch Screen
- Lung Mechanics
- Heated Expiratory Channel
- Spontaneous Patient Data with Graphs and Loops



T-Bird Legacy / Vs02 / AVS

- Transport / ICU / Home
- Volume Control, Pressure Control, SIMV, PS, CPAP & BiPAP
- Turbine Driven (No Compressor Required)
- Monitoring on Expiratory Parameters
- Internal Battery Backup
- F102 (Vs02 / AVS) Available



PB Make LP10

- Transport Category Ventilator
- Volume Control Ventilator
- CMV & SIMV Mode
- Pressure Limiting Ventilator
- Internal Battery Backup



Refurbished CPAP

- Auto Cpap And Manual Both Type Of
- Machines Are Available
- Optional Humidifier Available
- Full Leakage Compensation
- Resptronics And Resmed



Refurbished Bi-Pap

- Bi-level Spontaneous And Back Up Rate
- Both Types Of Machines Available
- Full Leakage Compensation
- Automatic Trigger Sensitivity
- Display For Expiratory Pressures
- Ramp Facility With Time And Pressure Adjustment



Accessories





MEDICAL EQUIPMENTS



A Division of Nira Group

506, B/3, Near AIIMS Hospital Tatibandh, Raipur (Chhattisgarh) - 492099

H-10 / H-100

- LCD Display
- Small & Lightweight Design
- SpO2, PR, Bar Graph and Plethysograph Display
- Easy to use, up to 300 hrs SpO2 data
- Std. 4 "AA" Batteries Power Supply & Rechargeable Battery
- Audible & Visual Alarms



M3A Vital Signs Monitors

- 3.5" high resolution TFT LCD & clear LED numeric Display
- SpO2 / SpO2 & NIBP parameters Measurements
- Adaptability to Adult, Pediatric & neonatal Patients
- Long Battery Period up to 21 hr. While generating less noise
- High Capability data storage & Nurse call
- Smart data Management Software for Network



SE-3 Three Channel ECG

- Portable & user friendly Design
- Foldable LCD Screen
- One Button Operation
- 12 Leads Simultaneous Acquisition
- Built In High Resolution Thermal Printer
- Rechargeable Li-ion Battery for Continuous 3 hrs Use
- 12 Leads Grid Printing on Normal Roll or A4 paper by External Printer
- Large Data Storage up to 144 ECGs & can be Magage by USB flash Disk
- Automatic ECG measurements & interpretation tested by CSE/AMAMIT database



SE-3 Twelve Channel ECG

- Portable & user friendly Design
- 5.7" LCD Screen with blue back light
- One Button Operation
- Alphanumeric Keyboard
- Rechargeable Li-ion Battery for Continuous 300 copy ECG print.
- 12 Leads Simultaneous Acquisition
- Built In High Resolution Thermal Printer & A4 paper by External USB Printing
- Large Data Storage up to 100 ECGs & can be Magage by USB flash Disk
- Automatic ECG measurements & interpretation tested by CSE/AMAMIT database



CM-3600 Portable Patient Monitor

- 8.4" Color TFT Display
- 3/5 Parameters Monitor
- Built in rechargeable battery up to 3 hrs backup
- Built in drug dosage calculation & Oxy-CRG mode
- Adaptability to Adult, Pediatric & Neonatal Patients
- 7 waveforms maximum with audible & visual alarms



Patient Monitor (C.E. Approved)

- 12.1" high resolution color TFT display with max 11 waveforms
- 3/5 lead ECG, NIBP, SpO2, RESP, Temp, PR
- Pacemaker detection, electrocautery interference proof
- Defibrillation protection and defibrillation synchronization
- Usb and serial Ports enables abundant future upgrades
- Oxy CRG Available to judge the respiration for neonates
- Large font display & Data Storage
- Arrhythmia and S-T segment analysis
- Rechargeable Li-ion (4 hours max.)
- Optional : 12 - lead ECG, Nellcor - SpO2, Omron NIBP 2-IBP / 4-IBP.
- Mainstream / sidestream C2
- 2-USB ports and SD card slot ensure data storage / function upgrades



Fetal Doppler

- Easy - to use
- Compact Design
- Accurate fetal heart detection
- Interchangeable Water Proof Probe
- Large Backlight LCD Display



Fetal Monitor F-6

- FHR, TOCO, Event Mark, AFM
- 10.1" High resolution color folding Lcd Screen
- Optional Dual FHR Monitoring
- 24 hours waveforms Playback
- 150/82 mm wide Paper printing
- 12.3 cm / min real time printing speeds
- Quick printing for stored waveform
- Lithium battery for 4 hours continuous working
- Probe rack and wall mounting rack
- Handle for easy carry
- Multi-crystal pulse wave doppler waterproof transducer for fhr detection



Syringe Infusion Pump

- Front loading of syringe
- Data log event : 500 data log events in real time
- Dynamic pressure system (DPS)
- Auto detection of syringe
- Works with standard and wide used brands of syringe
- Option to add customized and store brands up to 20 no
- Syringe capacity : 5 to 80 ml
- Flow rate range : up to 2000 ml / hr.
- Drug library : 64 drugs (user editable)
- Various mode of infusion
- Eg. Disconan, BD precise, plastipak, BBRAUN perfusor, romeo etc.



Drop Infusion Pump

- Keypad Locking
- Free flow protection
- Simplified selection
- Lighted status indicator
- Customized IV sets library
- Comprehensive & fine display
- 7-Segments Led & 4 Lensear LCD Display
- Option to select flow rate unit in terms of drop / min or ml / hr. & volume unit in terms of drops or ml



Fisher & Paykel Bubble CPAP

- Non invasive method of ventilation
- MR 850 Humidifier for optimal heat 37o
- Single heated birthing circuit
- Limits Pressure to 17 cm H2O at 8L / min
- Provides constant CPAP
- Easy to maintain & use
- Cose system for patient safety



Fisher & Paykel Airvo2

- Non invasive method of ventilation
- It provides high flow with air oxygen mixtures
- It has a built - in humidifier, Blender & air source
- It flows from 2 to 60 liters per min with nasal interface
- It has built - in disinfectant kit
- Dual spiral heated birthing tubing



Radiant Heat warmer

- Micro controller base skin / manual mode
- Provide thermal regulation with highest degree of accuracy
- User friendly feather touch control
- Auto calibration for temperature probe
- Clear acrylic collapsible side panels
- Available with resuscitation / fixed cradle / trolley type models



Photo Therapy Unit

Available Models :-

- LED / CFL / TUBE base double surface photo therapy
- LED / CFL / TUBE photo therapy stand
- LED / CFL / TUBE base under surface



Refurbished defibrillator

- Monophasic / Biphasic Defibrillator with recorder
- 3 Lead ECG Monitoring
- Reusable Adult & Pediatric paddles
- Adjustable alarm limits for HR
- manual & synchronize mode



Anesthesia Workstation

- Mild steel body with powder coated for high corrosion resistance
- S.S. Table top & Monitor stand on eye level
- Goldman Hiothane vaporizer



Anesthesia Workstation

- 3 GAS anesthesia delivery system
- Oxygen (Dual) Nitrous (Dual) air
- Highest level of delivery with min flow
- Support wide range of vaporization
- Top ray for mounting patient monitor
- Touch screen ventilator with volume, Pressure control and support with PEEP



Electro Surgical Unit's

- Models Available
- Micro-controller based mono - polar, vessel seilers
- Argon plasma unit



Surgical Operating Table

- Motorized remote controlled table
- C-ARM Compatible Table
- Basic hydraulic Table
- All ortho Attachments



LED Surgical Operating Light

- Shadow less daylight with color change
- Significantly no heat - up in operating field
- No tissue dry-up during extended operation
- Excellent depth illumination

